



Reagenssats, hårdhet (total), HA-71A testsats

Produktnr:

2437800

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For determination of Total Hardness by the Drop Count EDTA Titration Method. Replacement reagent set for HA-71A Test Kit. Set contains Buffer Solution (Hardness 1), Hardness 2 and Hardness 3 Test Solutions for approximately 100 tests.

Specifikationer

Mätområde:	1 - 20 mg/L CaCO ₃
Metodnamn:	EDTA Titration
Parameter:	Hardness (Total)
Specialfunktioner:	Reagent set, HA-71A test kit