



Reagent set, LR molybdenum, 5 mL

Produktnr: 2439500

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



For determination of Low Range Molybdate Molybdenum by the ternary complex method, using visual test kits. Replacement reagent set for MO-LR Test Kit (2359300). Range: 0.25 - 3 mg/L Mo. Sample Size: 5 mL. Set contains: Molybdenum 1 Reagent Powder Pillows and Molybdenum 2 Reagent Solution for approximately 100 tests.

Specifikationer

Mätområde:	5 mL
Metodnamn:	Ternary Complex