



Reagenssats, sulfit, SU-5 sats

Produktnr:

2440000

SEK Pris:
Tillgänglig

Kontakta oss

Farligt



Reagent set for determination of Sulfite by iodometric drop count titration. Replacement reagent set for SU-5 Test Kit (148002). Range: Up to 200 mg/L SO₃. Set includes Sulfamic Acid Powder Pillows, Sulfite 1 Reagent Powder Pillows and Sulfite 3 Solution for approximately 100 tests.

Specifikationer

Specialfunktioner:

Sulfite Reagent set, SU-5 Kit : Iodometric Method