



Hårdhet (kalcium och total sekventiell) HR-reagenset, digital titrator

Produktnr:

2448700

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For sequential determination of Calcium and Total Hardness by EDTA Digital Titrator titration. Hach method 8329. Range: 100 - 4000 mg/L CaCO₃. Set contains Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), Sulfuric Acid (5.25 N) and 0.800 M EDTA Digital Titrator Cartridge. Approximately 40 - 100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	100 - 4000 mg/L CaCO ₃
Metodnamn:	EDTA Titration
Metodnummer:	8329
Parameter:	Hardness (Calcium and Total Sequential)
Plattform:	Digital Titrator
Specialfunktioner:	Reagent set, HR, Digital Titrator