



Reagenssats, sulfit, titrering med burett

Produktnr:

2459800

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



For determination of Sulfite by iodometric buret titration. Hach Method 8071 Range: 0 to 500+ mg/L SO_3^{2-} . Reagent set includes: Dissolved Oxygen 3 Reagent Powder Pillows, Potassium Iodide-Iodate Standard Solution and Starch Indicator Solution for approximately 100 tests.

Specifikationer

Antal tester :	100
Mätområde:	0 - > 500 mg/L SO_3^{2-}
Metodnummer:	8071
Specialfunktioner:	Sulfite Reagent set, Buret Titration