



## Filter holder, easy pressure syringe, PALL

<b>Produktnr:</b>	<b>246800</b>
SEK Pris:	Kontakta oss
Tillgänglig	

Pack of 6 filter holders. Accepts 25 mm filter discs for fast filtration of 1-10 mL of sample.

**Fits most standard syringes with Luer-type fittings**

**Made of POM plastic with type-316 stainless steel membrane support screen**

**Autoclavable**

**Seal plugs included**

**Single filter PN 246801.**

---

### Specifikationer

Diameter, yttre:	25 mm
Enhet:	PK 6