



Be Right™



Flask, Erlenmeyer, wide mouth, 1000 mL, Subdivisions 250 mL

Produktnr: 2489453

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included.

Specifikationer

Enhet:	1 each
Kapacitet :	1000 mL
Material:	Glass Wide Mouth