



Sample Cell: 1" Square Glass 10 mL matched set of 8

Produktnr:	2495408
SEK Pris:	Kontakta oss
Tillgänglig	

Comes as a matched set of 8 cells.

Use in Combination with Hach Spectrophotometers

Hach spectrophotometers, chemistries, and methods are designed, developed, and manufactured to optimize your results.

Specifikationer

Fotnot :	1 pair of cells is included with purchase of DR 5000 UV-Vis, DR 3800, DR 2800, or DR 2700 Spectrophotometer
Volym/förpackningsstorlek:	8 /pk