



## Test kit, corrosion management, model CM-1

Produktnr: 2495700

### UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

To reduce lead and copper levels in distribution, service and domestic water lines, suppliers must stabilize aggressive water by adjusting pH, temperature, and hardness. Analysts can verify the success of their stabilization by screening samples with this rugged, portable test kit.

**Uses drop count titration and color disc methods**

**Maintain corrosion control of aggressive distribution system waters to meet Lead and Copper Rule Requirements**

**Establish baseline conditions to monitor treatment or disinfection strategy changes**

**Calculate Langlier and Agressive Indexes**

---

### Specifikationer

Antal tester :	100
Applikation:	Corrosion Management
Plattform:	Color Disc