



Parafilm M Laboratory Film

| | |
|-----------------------|---------------|
| Produktnr: | 251764 |
| SEK Pris: | Kontakta oss |
| Skickas inom 2 veckor | |

Seals labware quickly, effectively. Film prevents sample evaporation and contamination; protects anhydrous materials. Can be folded repeatedly and twisted around sharp surfaces without rupturing. Stretched lengthwise at 21 °C (71 °F), Parafilm M increases approximately 200%. Becomes more flexible, softer, stronger as temperature increases. Resistant to alcohol, medium concentration acid, and alkalies. Width: 10 cm (4"); 38 m (125') roll.

Specifikationer

| | |
|----------------------------|--------|
| Volym/förpackningsstorlek: | 1 each |
|----------------------------|--------|