



Filter, glass fiber, diameter 47 mm, preweighed standard methods TSS procedure, 100/pk

Produktnr:	2546100
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Specifikationer

Diameter:	47 mm
Enhet:	100/pk
Material:	Glass Microfiber
Porstorlek:	1.5 µm
Volym/förpackningsstorlek:	100 /pk