



Sample container, 45 mL

Produktnr:	2552114
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Plastic, leak-proof. Corning Snap Seal[™] offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

Specifikationer

Diameter, yttre:	30 mm (approx)
Kapacitet :	45 mL
Volym/förpackningsstorlek:	600 /pk