



Sample container, 120 mL, pk/250

Produktnr:	2552212
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Plastic, leak-proof. Corning Snap Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

Specifikationer

Diameter, yttre:	45 mm (approx)
Kapacitet :	120 mL - 250/pk
Volym/förpackningsstorlek:	250 /pk