



Test kit, chloride, model CDS-DT

Produktnr:

2580600

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikationer

Antal tester :	100
Kapslingstyp:	F
Mätområde:	10 - 10,000 mg/L Cl ⁻
Metodnamn:	Digital Titrator/ Silver nitrate
Parameter:	Chloride - As Cl ⁻ , 10-10,000 mg/L Cl ⁻
Plattform:	Digital Titrator
Smallest increment:	10
Test-kit modell:	CDS-DT

Obligatoriska tillbehör

- Pipett, variabel, volym 1,0 - 5,0 mL (Item BBP065)