



Be Right™



Filter assembly, portable, complete (millipore microfil)

Produktnr: 2586200

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Föreslagna ersättningar:

- 2586201

For membrane filtration features unique membrane lifting device. Depressing the lever breaks the vacuum and lifts the membrane above the 47 mm membrane support for easy retrieval with sterilized forceps. Designed for use with Microfil funnels (see Filtration Funnels), creating an efficient, low-maintenance system for membrane filtration in the lab or field. Rugged metal construction with vacuum side-arm at the base. Works with any vacuum source.

Specifikationer

Volym/förpackningsstorlek:

1