



Reagenssats, hårdhet (ULR), bulklösning

Produktnr:

2603101

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Reagent set ultra low range hardness determination by the Chlorophosphonazo method. Hach Method 8374. Range: 4 - 1000 ug/L Ca & Mg as CaCO₃. Chlorophosphonazo reagent solution is packaged in a 500-mL bottle that is ideal for use with Hach's Pour-Thru Cell. Approximately 250 tests

Convenience

Rapid Liquid methods utilize liquid reagents, bottletop dispensers and pour-thru cells to provide simplified testing.

Performance

Since the pour-thru cell provides a stationary optical system, Rapid Liquid methods typically provide better precision than methods that require movement of sample cells.

Efficiency

Rapid Liquid methods are ideal for high-volume testing.

Specifikationer

Antal tester :	250
Instrument:	All except DR 900, PC II
Mätområde:	4 - 1000 µg/L CaCO ₃
Metodnamn:	Chlorophosphonazo
Metodnummer:	8374
Parameter:	Hardness
Plattform:	Lösning
Specialfunktioner:	Reagent set, ULR, Bulk