



Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

Produktnr:	2634001
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. (**Please Note:** These sample vials are not autoclaveable)

Specifikationer

Diameter:	43mm
Kapacitet :	fill line at 100 mL (+/-2.5 mL)
Längd:	93mm
Specialfunktioner:	Sample Vial, with Thiosulfate
Unit:	200 /pk
Volym/förpackningsstorlek:	200 /pk