



Buret, Class A, Glass, 25 mL

Produktnr:	2636540
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Glass, funnel-fill buret. Certified Class A (ASTM Specification E287) with permanent serial number and certificate for NIST traceability. PTFE stopcock does not require grease. Dust cap included. 10 mL buret meets requirements for potentiometric titrations (ASTM D664).

Specifikationer

Kapacitet :	25 mL
Material:	Class A Glass
Specialfunktioner:	Buret