



Be Right™



Gasolflaska, 1200 mL

Produktnr:	2662200
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

For use in precipitating elements from gas. Used in the Hach Cold Vapor Mercury Apparatus Set (Prod. No. 2674400) to perform Cold Vapor Mercury Concentration Analysis—Hach Method 10065.

Specifikationer

Body Material:	Glass
Kapacitet :	1200 mL
Material:	Borosilicate Glass
Package size UOM:	1pk
Volym/förpackningsstorlek:	1 each