



Be Right™



## Test kit, hypochlorite, model CN-HR

Produktnr:

2687200

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Hypochlorite Kit, Model CN-HR - Drop Count

Utilizing the simple, economical drop count titration method, this kit is used to determine high range hypochlorite solution concentrations (bleach solutions). Range: 5-15%, contains 100 tests.

**Simple, economical titration method**

**Kit contains all required reagents, solutions, and glassware**

### Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	50 - 150 g/L Cl <sub>2</sub> (5-15%)
Metodnamn:	Drop count titration/
Parameter:	Thiosulfate
Plattform:	Hypochlorite, high range - (Bleach) As Cl <sub>2</sub>
Smallest increment:	Drop Count Titration
0.5%	
Test-kit modell:	CN-HR