



Vacuum filter microfunnel, pk/50.

Produktnr:	2831500
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Feature a completely assembled filtration unit: Funnel, filter, nutrient pad and petri dish. Filter: GN-6 Metrical with 0.45 µm pore size and 47 mm diameter; white, gridded surface. Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Nutrient Pad: Cellulose. Fits 2831500, 2846100 and 2846150 manifolds. Sterile and individually bagged.

Specifikationer

Diameter:	47 mm
Material:	Mixed-Cellulose Ester
Porstorlek:	0.45 µm
Volym/förpackningsstorlek:	50