



Membrane filters millipore S-pak type HA, 47 mm, 0.45 um pore size, pk/600

Produktnr:	2936500
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 um. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Specifikationer

Diameter:	47 mm mm
Material:	Mixed-Cellulose Ester
Porstorlek:	0.45 µm
Volym/förpackningsstorlek:	600