



Fast volym, mikropipett, 1000 µL

Produktnr:	2946010
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and color-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specifikationer

Kapacitet :	1000 µL
Noggrannhet:	0,5 %
Volym/förpackningsstorlek:	1 each