



MultiTest kit case (295 x 330 x 576 mm), blue polypropylene

Produktnr:	4660800
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Durable blue polypropylene carrying case allows easy transport of items for field testing.

Recessed handle in lid provides easy stacking

Contents are carried in upright position

Includes recessed flat-panel areas for labeling

Specifikationer

Dimension (H x B x D):	295 mm x 330 mm x 576 mm
Höjd :	11 in
Kapslingstyp:	J