



## OriFlo™ Filter apparatus

<b>Produktnr:</b>	<b>4966000</b>
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Designed for rapid, qualitative filtration of particulate matter from aqueous samples using a 25 mm glass fiber filter or other medium with porosity greater than or equal to 1 micron.

**Used in Hach Method 10014 for Ultra Low Range Chlorine in Wastewater**

**Apparatus included in the Ultra Low Range Chlorine Apparatus Set, PN 2595600**

**Not designed for use as a microbiological filtration device using membrane filters with a porosity of 0.45 micron or less**

**Cannot be autoclaved or sterilized using organic solvents**