



## Beaker, Griffin, Low Form, Glass, 600 mL

<b>Produktnr:</b>	<b>50052</b>
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specifikationer

Enhet:	1 each
Kapacitet :	600 mL
Material:	Glass, Griffin, Low Form