



Beaker, griffin, low form, glass, 100 mL, 12/pk

Produktnr:	50072H
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Specifikationer

Enhet:	12 /pk
Kapacitet :	100 mL
Material:	Glass, Griffin, Low Form