



## Cylinder, graduated, 100 mL +- 0.6 mL, 1.0 mL divisions (white markings)

|                       |              |
|-----------------------|--------------|
| <b>Produktnr:</b>     | <b>50842</b> |
| SEK Pris:             | Kontakta oss |
| Skickas inom 2 veckor |              |

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes  $\geq 25$  mL.

---

### Specifikationer

|                            |               |
|----------------------------|---------------|
| Certifieringar :           | Class B       |
| Höjd :                     | 255 mm        |
| Kapacitet :                | 100 mL        |
| Material:                  | Glass         |
| Mätmetod:                  | To Deliver    |
| Typ:                       | Single Metric |
| Volym/förpackningsstorlek: | 1 each        |