



Container, glass, 10 L (2.5 gal), with PTFE-lined cap

Produktnr: 6559
SEK Pris: Kontakta oss
Inget leveransdatum anges

2.5 gal glass container with cap

Specifikationer

Vad ingår i leveransen:

Bottle(6358)

Cap (6355)

PTFE Liner (6378)

Add Required parts:

For Compact or Standard 12 Bottle Portable Base

Container Support (1502)

For Refrigerated and All Weather Samplers

Container Support (1502)

Extension Tube (3527)

Vad finns i förpackningen

Bottle(6358) Cap (6355) PTFE Liner (6378)Add Required parts:For Compact or Standard 12 Bottle Portable BaseContainer Support (1502)For Refrigerated and All Weather SamplersContainer Support (1502)Extension Tube (3527)