



Pensel, burett, 25 mL

Produktnr:	68500
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Ideally shaped for quick and thorough cleaning of burets and other similarly shaped vessels. Unbleached white nylon bristles are supported by twisted galvanized steel wire. Brush portion measures 8 cm (3.1") in length. Fan style tip. 1.3cm (0.5") Dia. 76cm (30") L

Specifikationer

Volym/förpackningsstorlek:	1 each
----------------------------	--------