



Radiometer Analytical XR110 referenselektrod (röd stav, jordanslutning, skruvlock)

Produktnr: B20B112
SEK Pris: Kontakta oss
Inget leveransdatum anges

The XR110 Reference Electrode has a glass body, Red Rod Ag/AgCl reference technology system, refillable electrolyte (KCl 3M), ground joint, and porous pin junction. The XR110 reference probe is recommended for Voltammetry applications and Educational purposes. The XR110 is a screw cap electrode and requires an external connection cable for connection to a meter. (Replace Obsolete #B20B110)

Specifikationer

Dimension (D x L):	8 mm x 120 mm
Elektrolyt:	KCl 3M
Garanti :	12 månader
Givar koppling:	Porous plug
Givartyp:	Laboratory reference electrode
Kontakt:	Screw cap S7
Material:	Glass
Temperaturområde :	-10 to 100°C
Vad ingår i leveransen:	Radiometer Analytical XR110 referenselektrod, 30 mL KCl 3M påfyllningslösning och bruksanvisning.

Vad finns i förpackningen

Radiometer Analytical XR110 referenselektrod, 30 mL KCl 3M påfyllningslösning och bruksanvisning.

Obligatoriska tillbehör

- CL111 Elektrod kabel FX/FIL/1m/B (Item A94L111)
- CL114-kabel, 1m, för skruvlockselektrod FX/S7/coax, BNC-kontakt (Item A94L114)