



Sedimenter for AISE, NISE and AN-ISE sc

Produktnr:	LZY771
SEK Pris:	Kontakta oss
Skickas inom 1 vecka	

The sedimenter can be used as a bypass solution. The sample can be let through the sedimenter without prefiltering.

Specifikationer

Handboksspråk:	Non verbal instruction sheet
Instrument:	suitable for: all AISE sensors; all NISE sensors; all AN-ISE sensors.
Längd:	Intallation plate for wall mount: 800 mm x 495 mm
Material:	PVC
Monteringsform:	Wall installation (indoor only).
Vikt :	Sedimenter incl, wall plate only: 14,7 kg