



Be Right™



Filterprobe sc, 10 m heated hose, 230 V

Produktnr: LXV429.99.02000
SEK Pris: Kontakta oss
Tillgänglig

In-situ membrane filtration. The immersed filter probe takes samples continuously. Air bubbles prevent coatings from forming on the filter module.

Specifikationer

Dimension (H x B x D):	384 mm x 330 mm x 243 mm
Flöde:	≥ 5 mL/min filtrate flow rate, filter probe pumps 4 out of 5 minutes
Kapsling:	Plastic enclosure, PPE, flammability class in accordance with UL 94
Lagringsförhållanden:	-20 °C - 60 °C
Neddopningsdjup:	måx. 3.0 m
Omgivningstemperatur:	-20 - 45 °C sample tubing
pH-område:	4 - 9
Provtemperatur:	4 - 40°C (39 - 104°F)
Säkringsskydd:	Via analysis instrument, via controller
Slanglängd:	10 m heated
Specialkommentarer:	Components: filter module wih pump; once the filter modules have become wet, they must be stored wet
Vikt :	8kg hose (10 m)