



## Total Alkalinity Test Strips, 0-240 mg/L, 50 tests

Produktnr:	2744850
SEK Pris:	Kontakta oss
Tillgänglig	

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to drinking water, food, and beverage processors; boiler operators; and aquatic biologists. Contains 50 tests.

**Easy to use, disposable, and inexpensive**

**A great way to obtain quick, quantitative answers in the field or in the lab**

---

### Specifikationer

Antal tester :	50
Enhet:	50 tests
Mätområde:	0 - 240 mg/L
Minsta ökningssteg:	Steps: 0, 40, 80, 120, 180, 240 mg/L (ppm)
Parameter:	Total Alkalinity - Phenolphthalein
Plattform:	Teststickor
Test-kit modell:	Test Strips