



Glass Fiber Filter, Diameter 70 mm, Preweighed for Standard Methods TSS procedure, 100/pk

Produktnr:	2546135
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare

Specifikationer

Diameter:	70 mm
Enhet:	100/pk
Material:	Glass Microfiber
Porstorlek:	1.5 µm
Volym/förpackningsstorlek:	100 /pk