



Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

Produktnr:	2096250
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Specifikationer

Diameter:	25 mm
Dimension (D x L):	25 mm x 150 mm
Kapacitet :	Capacity to Neck 50 mL
Längd:	150 mm
Specialfunktioner:	Tube, 25 x 150 mm
Volym/förpackningsstorlek:	10 /pk