



Be Right™



Manifold, filter funnel, aluminum

Produktnr: 2486100
SEK Pris: Kontakta oss
Ring för leveransdatum

Feature 3-place filtration with independent vacuum control. **Adapters are not included**, but are ordered separately from PALL. (Adapters shown in the photo not included. These are elongated standard adapters)

Accepts 25 or 47 mm filter funnels for processing multiple samples simultaneously. Designed for efficient analysis of microbiological or particulate contamination. Aluminum: Anodized aluminum body; chrome-plated brass valves; stainless steel drain plugs. Autoclavable. Dimensions (H x W x D) 3.9" x 17.2" x 6.0"

Specifikationer

Volym/förpackningsstorlek: 1 each