



Bag, Whirl-Pak* 5-1/2X9, pk/10

Produktnr:	1437297
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection.

Specifikationer

Kapacitet :	720 mL
Volym/förpackningsstorlek:	10