



Be Right™



Buffer, magnesium chloride/potassium phosphate, 99 mL, pk/25

Produktnr: 1430598
SEK Pris: Kontakta oss
Ring för leveransdatum

Magnesium Chloride/Potassium Phosphate Buffer, 99 mL. Buffered Dilution Water (required for testing non-potable water). pk/25.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Specifikationer

Volym/förpackningsstorlek: /pk