



**Be Right™**



## **MEL/P-A safe drinking water lab**

**Produktnr: 2569600**

**UTGÅNGEN ARTIKEL**

Den här artikeln är inte längre tillgänglig.

### **MEL/P-A Safe Drinking Water Lab**

Economical portable lab monitors key water quality parameters plus Total Coliforms and E. Coli (Presence-Absence Method). MEL P/A Safe Drinking Water Lab includes portable incubator, reagent sets, portable long-wave UV lamp, manual and apparatus.

#### **Consumables only**

**Choose media and order seperately**

**Nitrate and Chlorine determined using visual color discs**

**Portable incubator used to determine total and fecal coliforms directly in the field**

---

### **Specifikationer**

Känslighet:

CFU/mL