



MEL/850 Portable,

Produktnr: 2688800

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

MEL/850 Portable Water Lab

MEL/850 Portable Water Laboratory includes: DR/850 Colorimeter, Pocket Pal TDS and pH Testers, Portable Incubator, Portable UV Lamp, media and apparatus for testing E. coli and total coliforms, plus 10 common inorganic water quality tests.

DR/850 Colorimeter with reagents to monitor 8 parameters

50 additional parameters available on DR/850 Colorimeter; order reagents seperately

Consumables and P/A media included

Portable incubator used to determine total and fecal coliforms directly in the field

Specifikationer

Känslighet: 1 CFU/mL