



ComLabo Easy Datorprogram för Sension+ GLP, CD, kabel, USB-adapter

Produktnr:

LZW8997.99

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

LabCom Easy software, to obtain data of Sension+ GLP instruments from a PC. Kit including CD and cable with 9-contact connector and USB adapter.

- On line graph view of measured value
- Results table personalisation, the user selects the displayed information
- Storage of calibration and measuring results in a data base where the data can be viewed or printed at any moment
- Print report of individual sample or results summary for several samples

Data export to Excel

Convenient Data Management.

Labcom easy allows the acquisition and the management of the obtained results from Sension+ Laboratory instruments.