




## SAC850 Sample changer, Standard version with cables

**Produktnr:**

**R41T043**

### UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Basic sample changer for automation of small sample series, single-radius sample trays for up to 24 samples.

**Single-radius sample trays for up to 24 samples**

**Parking position on sample tray**

**Dip rinsing/electrode conditioning performed on 1 or 2 dedicated positions on sample tray**

---

### Specifikationer

|                        |                          |
|------------------------|--------------------------|
| Dimension (H x B x D): | 711 mm x 580 mm x 686 mm |
| Längd:                 | 686 mm                   |
| Typ:                   | Provväxlare              |