



TEST KIT, UNIV DIGITAL TITRATOR

Produktnr:

2270900

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Accurate, convenient titrations without the bulk, fragility, or waste of a conventional buret.

Universal Digital Titrator Kit with Digital Titrator, manual, 100-ml graduated cylinder, 125 and 250 ml Erlenmeyer flasks, demineralizer bottle, clippers, and rugged carrying case.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40 - 100 typical titrations per cartridge.

Specifikationer

Leveransomfattning:	Digital Titrator with accessories in carrying case
Noggrannhet:	± 1 % (Uncertainty of readings is 1 digit. Most samples require more than 100 digits)
Tillbehör som ingår:	Ja
Vad ingår i leveransen:	Includes Digital Titrator, manual, 10-mL graduated cylinder, 125 - and 250-mL Erlenmeyer flasks, demineralizer bottle and a carrying case.
Vikt :	0,132 kg

Vad finns i förpackningen

Includes Digital Titrator, manual, 10-mL graduated cylinder, 125 - and 250-mL Erlenmeyer flasks, demineralizer bottle and a carrying case.