



5500sc Underhållssats för ammoniak och monokloramin, 1 kanal

Produktnr: 9560501
SEK Pris: Kontakta oss
Inget leveransdatum anges

This comprehensive maintenance kit is designed to support the upkeep of your 5500sc Ammonia Monochloramine Analyzer, ensuring consistent performance and minimal downtime. It includes all necessary components for routine servicing and replacement.

Kit Contents:

Tubing & Fittings

Silicone, PVC, TPE, and Tygon tubing in various sizes (e.g., 1.4 m, 3.23 m, 1.04 m, 1.37 m, 0.37 m)

Multiple fittings including elbows, tees, reducers, barbs, adapters, and bushings (1–3 units each)

Filtration & Flow Control

Air filter (10 micron) – 1 unit

Liquid filter (40 micron, 100 PSI) – 1 unit

Flow restrictor – 1 unit

Check valve – 1 unit

Pump & Leak Detection

Liquid pump assembly with bracket – 1 unit

Leak detector circuit board – 1 unit

Colorimeter Components

Capillary tubing assemblies (red and black) – 3 units each

Stir bars (Teflon, 8x3mm) – 2 units

Cell caps – 2 units

Reagent Bottle Caps

Orange, blue, green, red, yellow, black – 1 unit each

Documentation

Tubing replacement instructions

Drain and vent tubing replacement instructions

Miscellaneous

Screws, thread tape, labeling materials, and packaging components for secure installation and maintenance

This kit provides all the essential components for maintaining the analyzer's fluid handling, measurement accuracy, and operational integrity.