



5500sc Ammoniak Monokloramin Standard 1 (0 ppm NH₃), 2 L

Produktnr:	25237000
SEK Pris:	Kontakta oss
Tillgänglig	

Standard 1 for 5500sc AMC Analyzer, 2 L

Specifikationer

Instrument:	5500sc
Lagringsförhållanden:	Keep container tightly closed in a dry and well-ventilated place. Keep in properly labeled containers.
Mätområde:	N/A
Parameter:	Ammonia/Monochloramine

Obligatoriska tillbehör

- Skyddshandskar storlek L, blå, nitril, puderfria, 50 styck (Item SM743L)