



Be Right™



## 5500sc Ammoniak monokloramin surt ytaktivt tvättmedel, 2 L

Produktnr:

25239000

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Cleaning solution for 5500sc AMC Analyzer, 2 L

---

### Specifikationer

Instrument:	5500sc
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Mätområde:	N/A
Parameter:	Ammonia/Monochloramine

---

### Obligatoriska tillbehör

- Skyddshandskar storlek L, blå, nitril, puderfria, 50 styck (Item SM743L)