



Reagenssats för klorid, 0,1-25,0 mg/L Cl⁻

Produktnr:

2319800

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Reagent set for the determination of Chloride, Mercuric Thiocyanate method.

Ej tillgänglig i Sverige pga. kvicksilverförbudet

Specifikationer

Antal tester :	100
Fotnot :	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Instrument:	All except DR 900
	DR/800 Series, PC II
Mätområde:	0 - 25,0 mg/L Cl ⁻
Metodnamn:	Mercuric Thiocyanate
Metodnummer:	8113
Parameter:	Chloride
Plattform:	Lösning