



Jodat-jodid, Digital Titrator-patron, 0,3998 N

Produktnr:

1496101

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Bruksfärdig titreringslösning för analys av sulfit, 4-400 mg/l SO₃.

Iodate-Iodide Standard Solution. For determination of sulfite by Digital Titrator iodometric titration. Digits = mg/L SO₃ in 20 mL sample.
Digits x 0.0005 = meq. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Mätområde:	0 - 800 mg/L SO ₃ ²⁻
Metodnamn:	Iodate-Iodide
Parameter:	Sulfite
Plattform:	Digital Titrator
Titrant:	0.3998 N KIO ₃ -KI