



Chlorine (free & total) colour disc test kit, model CN-66 (pk/100)

Produktnr:

223101

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Measuring range 0.1-3.4 mg/L Cl₂

Specifikationer

Antal tester :	50 Free, 50 Total
Kapslingstyp:	B
Mätområde:	Single Range
Metodnamn:	Color disc/DPD
Modell:	CN-66
Parameter:	Chlorine (Free and Total)
Plattform:	Color Disc
Smallest increment:	0.2
Standardmätområde:	0 - 3.4 mg/L Cl ₂
Test-kit modell:	CN-66 ²